

A.3500

A.3525

A.3510

A.3510.100

LUJ

gel lubrifiant lubricating jelly

Transparent - Hydrosoluble

Transparent - Watersoluble



CE
0123

CE

LUJ de Sapi Med est un gel transparent, neutre et hydrosoluble qui exerce une lubrification prolongée et continue. Il est compatible avec les muqueuses, les outils endoscopiques flexibles et les instruments chirurgicaux. La gamme propose plusieurs emballages et caractéristiques pour mieux répondre aux exigences du spécialiste et du patient.

LUJ from Sapi Med is transparent, water-soluble and sterile, and provides continuous and prolonged lubrication. It is compatible with the mucous membranes, flexible endoscopic instruments and surgical instruments. The line offers a wide range of packages and formats to satisfy the needs of both specialist and patient.

Ingrédients Ingredients:

Eau / water • Propylene Glycol • Carbomer • Phenoxyethanol • Imidazolidinyl Urea • Sodium Hydroxide • Methylparaben • Disodium Edta • Ethylparaben • Butylparaben • Propylparaben.



LUJ

gel lubrificante lubricating jelly

Transparent - Hydrosoluble
Transparent - Watersoluble

W.000118IGB REV 02 04/2008

REF A.3500



Distributeur de 500 ml en plastique recyclable.

Gel **non stérile** pour explorations rectales, gynécologiques et lubrification d'outils endoscopiques flexibles.

500ml dispenser in recyclable plastic.

Non sterile jelly for rectal and gynaecological examinations, and to lubricate flexible endoscopic instruments.

Emballage: distributeur de 500 ml

Package: 500ml dispensers

REF A.3525



Tube de 250 ml. Gel **non stérile** pour explorations rectales, gynécologiques et lubrification d'outils endoscopiques flexibles.

250ml tube. **Non sterile** jelly for rectal and gynaecological examinations, and to lubricate flexible endoscopic instruments.

Emballage: tube de 250 ml

Package: 250ml tube

REF A.3510



Sachet unidose de 10 ml. **Gel stérile** pour cathétérisme, explorations rectales, gynécologiques et lubrification d'outils endoscopiques flexibles.
Etui de 20 sachets.

10ml single use sachet **sterile** jelly for catheter insertion, rectal and gynaecological examinations, and to lubricate flexible endoscopic instruments. Box of 20 sachets.

Emballage: étui de 20 sachets

Package: box containing 20 sachets

Emballage pour les hôpitaux:
A.3510.100 : boîte de 100 sachets

Package for hospitals use only:
A.3510.100: box containing 100 sachets

Stérilisé par: rayons bêta

Sterilised: with beta radiation